

Adjust. Inject. Perfect.

Precision in every detail.

The new benchmark in premium needles for professionals.





"Now adjust depth with precision—right at your fingertips."

Achieve Precision in Every Procedure

Unlike conventional needles with fixed lengths, the rotational hub structure enables real-time depth adjustment from 0–4 mm.

By tailoring the depth to match varying skin thickness across treatment areas, practitioners can perform procedures with greater precision and stability—enhancing delivery efficiency as well as overall injection comfort and patient satisfaction.

Whether for fillers, botulinum toxin, or various aesthetic injectables, this adjustable-depth needle provides the ideal solution for any aesthetic procedure requiring precise depth control.

Product Line-up

0.						
Gauge	34G	33G	32G	31G	29G	27G
Color						
OD	0.178	0.203	0.229	0.254	0.324	0.400
ID	0.110	0.110	0.150	0.150	0.230	0.310
ISO	ETW	TW	UTW	ETW	TW	TW
Remark	X NEW	X NEW	X NEW			



Precise Depth Control

Personalized, Area-Specific Delivery



Enhanced Drug Efficacy

Accurate Delivery with Minimal Waste



Reduced Procedure Time

Quick Workflow with No Tip Replacement

Inquiries

TEL: 031-529-8182 | FAX: 031-529-8184 | WEB: WWW.WOORHI.COM E-MAIL: WOORHIMED@GMAIL.COM

Woorhi Medical Co., Ltd.

23, Palya Sandan-ro 3beon-gil,

Jinjeop-eup, Namyangju-si, Gyeonggi-do, Korea




Unique features of Johns Screw Needle

- The world's first needle length controllable type (adjustable range: 0~ 4 mm)
- Ideal for vertical injection in target dermis layer
- Korean and PCT patent application filed
- Mainly used for HA fillers and Meso products injection
- Suitable for all types of syringes (up to 10ml size)

Vertical

Vertical



Specification of Johns Screw Needle

whole lenth	Needle length	Adjust.range	Types	Raw materials	
	7mm	0~4mm	27G(Grey) 29G(Red) 31G(White)	Surgical stainless	
30mm			32G(Deep Green) 33G(Black) 34G(Orange)	Steel (needle) PP (Hub)	

Caution for use

- Do not use if the package is damaged.
- Discard the needles into a suitable container after use.
- Single use only. Do not reuse.
- The needle cannot be re-sterilized for any reason, even not been used.









27G x 7mm (adjustable range 0~4mm) 100 Needles



29G x 7mm (adjustable range 0~4mm) 100 Needles



31G x 7mm (adjustable range 0~4mm) 100 Needles



32G x 7mm (adjustable range 0~4mm) 100 Needles



33G x 7mm (adjustable range 0~4mm) 100 Needles



34G x 7mm (adjustable range 0~4mm) 100 Needles