




DQ

Beautiful Performance



KFDA  CE

How to Use



DermaQueen

Redefining Aesthetic & Anti-aging Injection



OLED Control Display

- Dosage
- Speed
- Depth



Automatic Negative Pressure Release System

Hose & Filter

The best way to inject easily & safely

Simple Setting
Easy Operation
Safe Treatment


Comfortable Grip
Through Ergonomic Design

Decrease Medicine Loss & Unnecessary Waste


Our special plunger holder is designed to decrease medicine loss by pulling it back after injection

The best way to inject medicine to target layers equally


The more accurate and equal, the more effective



Holds the skin for easier injection



Lifts up the skin to avoid injecting danger area



Equal (Balanced) medicine delivery into the right layer (Dermis)



Full face



Neck



Breast



The back of the hand


Multi-Application Platform

- Hyaluronic Acid
- Toxin
- PRP
- Placenta
- Growth Factor
- Vitamin etc

Handpiece Holder
having Battery
Charging Function



OLED Touch Screen
Controlling The Negative Pressure




Multi Needle

31G & 32G
5pin, 9pin, 16pin

- Fast and accurate treatment possible
- Minimize pain and drug loss in the treatment area
- Available in combination with regular syringes and DermaQueen(Injector)



**DermaQueen
(Injector)**



Syringes



**Package
(36ea / 1box)**


Specification



AP 7.4V 2S 2,000mAh HD 15C
310 x 205 x 50
590g



AC100-240V 50/70Hz
380 x 280 x 240
6.8kg(Portable)/10kg(Stand)



31G / 32G
5pin, 9pin, 16pin

Package

Portable : 10kg (450 x 300 x 480)
Stand : 16kg (520 x 370 x 920)

Manufacturer

Woorhi Medical Co.,Ltd.

Address: 23, Paryasandan-ro 3beon-gil,
Jinjeop-eup, Namyangju-si, Gyeonggi-do,
Republic of Korea

Tel: +82 31 529 8182, Fax: +82 31 529 8184

www.woorhi.com / www.dermaqueen.co.kr
woorhisutures@gmail.com



310 mm

850 mm

500 mm



1,055 mm

For the improvement of injection techniques and development of various treatments, DermaQueen is prepared to make constant efforts for new innovations.



www.dermaqueen.co.kr